STARTLE/PPI SYSTEM *For psychosis, schizophrenia and sensorimotor studies.*



b				UB Startle PP					
EXPERIMENT	PROTOCOL	RUI	VIE	W CONTR	OL CALIBR	ATE			
Export Export Open Output Selected All Folder CSV/Excel File	ut Deleted Delete Selected All Delete	Selected Experiment (ected ent (brief)					
Experiments	Experiment Info								
PPI_20211217_164755	All cubicles totals Pulse Alone	AVG PPI% 9 54,4	td.Dev. Std.E 20.93 9.	Pulse Alone		-			
PPI_20211217_162929	Pulse Alone Prepulse-Pulse	54.4	20.93 9.	Prenulse-Pulse		· · · · · · · · · · · · · · · · · · ·			
PPI_20211217_160657	PrePulse Alone	0	0	0 PrePulse Alone - No Stimulus					
PPI_20211217_152834	No Stimulus	0	0	0 No Stimulus	20	40	60	80	100
PPI_20211217_120606	PERCENTAGE OF PPI	3.1		U	20	40	60	80	100
PPI_20211217_110612	# Cubicle +	Timer	StepName 🖵	Stimulus 👻	Session 🔎	Protocol +	Experiment +	Treatmennt 👻	Stage
PPI_20211217_101651	6	182.979	Phase B	PrePulse Alone	PPI_20211217_16475	New BIO Protocol wit	Test experiment	Treated	Habitu 个
PPI_20211217_092817	7	193.855	Phase B	Pulse Alone	PPI_20211217_16475	New BIO Protocol wit	Test experiment	Treated	Habitu
PPI 20211217 092239	8		Phase B	No Stimulus		New BIO Protocol wit		Treated	Habitu
PPI 20211216 151512	9		Phase B	Prepulse-Pulse		New BIO Protocol wit		Treated	Habitu
PPI 20211216 150548	10		Phase B Phase B	Prepulse-Pulse PrePulse Alone		New BIO Protocol wit New BIO Protocol wit		Treated	Habitu Habitu
PPI_20211216_150324	12		Phase B	Prepulse-Pulse		New BIO Protocol wit		Treated	Habitu
PPI 20211216 145611	13		Phase B	PrePulse Alone		New BIO Protocol wit		Treated	Habitu 🗸
PPI 20211216 144023	40 record(s) found on 1.229 Elaosed time 00:00.48								
PPI 20211216 143941									
PPI 20211216 135754			START IPI		ld L	Da	GND		
PPI 20211215 163609	Stimuli [V]	_	STAR		START	Children of Stringson	DET		<u> </u>
PPI 20211215 145151						STIL			
PPI_20211215_145151	0.18 0.1	0.22	0.24	0.26 0.28	0.3 0.32	0.34 0.36	0.38 0.4	0.42 0.44	0.46
DDI 20211215_134420					Time [s]				





Pre-Pulse Inhibition of Startle Reflex (PPI)

Complete set up, ready to go.



Background

- The Startle/PPI system is a valuable tool to investigate complex human disorders, such as schizophrenia, which shows a dysfunction in the sensorimotor mechanism, as other forms of phycosis.
- Sensorimotor gating refers to the ability of a sensory event to suppress or decrease a motor response, in this case, a startle response.
- · The pulse alone startle reflex provides an index of

the animal reaction to the startle pulse in terms of maximal peak, the pre-pulse alone startle reflex provides an index of the animal reaction to the pre-pulse startle i and the %PPI indicates the inhibition of the startle reflex, due to the pre-pulse (in healthy animals).

• The %PPI is calculated as the decrease of the startle response after the pre-pulse is delivered, divided by the startle amplitude and is calculated for each of the pulse intensities).

Typical device applications

Several clinical studies have shown that a number of human disorders have impaired PPI including: schizophrenia, Huntington's disease, fragile X syndrome, and autism. Utilization of Startle PPI tests in rodents has a translational value to assess the effects of a number of treatment modalities such as putative anti-psychotics, and to explore genetic and neurobiological mechanisms underlying behaviors of relevance to psychosis. PPI deficits have been shown also in obsessive compulsive disorder, Tourette'ssyndrome and post traumatic stress disorder (Braff et al. 2001).

Detrimental effects on motor and sensory gating have also been shown to be associated to synthetic cannabinoids (Corli et al. 2023), to substances like ketamine (MArti et al. 2021) and the psychoactive effect of LSD impairs the correct brain processing (PPI) (Tirri et al. 2022).

Product Description

The Ugo Basile system is fully automated both for the preparation of the procedure and for the analysis of data. Main features include:

- Automatic detection of Startle Reflex in mice and rats.
- A typical system combines the following elements:
 - Dedicated Software;
 - Up to 4 Isolation Cubicles that include an infrared (I.R.) light, a loudspeaker and a silent fan (under 50 dB), all conveniently positioned inside the

sound attenuating box;

- Stimulating/Recording (S/R) Platform (includes the light, speaker and detection system using a force sensor).
- Different animal holders to match animal dimension.
- Flexible testing, easy to set up. Simply define trial number, sound and light stimulus and timing of experimental sequences (all fully randomizable):
 - Pulse; Pre-pulse; Inter-Pulse Interval; Inter-Stimulus Interval.

Procedures and Data Analysis

Pulses and pre-pulses consist of white noise and frequency sounds, lights or air-puffs.

The animal is restrained in a Plexiglas tube sitting on a force sensor. The dedicated software analyses data automatically, as startle response amplitude and %PPI.


The various procedures expose the animals to "Trials" of different nature:

- Pulse alone trial: only the strong stimulus is presented
- Pre-pulse alone trial: only the weak stimulus is presented
- Pre-pulse-pulse trial: the strong stimulus is preceded by the weak (pre-pulse) stimulus
- No-stim trials: no stimulus is delivered but the response of the animal is recorded to serve as a baseline

Experimental parameters can all be adjusted, linked, spaced in time, etc., to virtually build in a seamless and user-friendly manner virtually any type of startle and startle/PPI protocols.



Results and data are shown in a very intuitive way, as a summary data in the proprietary interface and as raw data (.csv), for those who wish to use their own analysis systems



Features	Benefits
User-friendly software included	Exceptional flexibility and intuitivity in the software protocol building
Configure your own Experimental Schedules directly on your PC	Maximum Flexibility, full event randomization
Multiple stimuli	The user can decide to use pre-pulses and pulses with white noise, frequency noise, light, air puff
Multiple cage system	High throughput with up to 4 cages running at the same time

Main references

G. Corli, M. Tirri, S. Bilel, R. Arfè, T. Coccini, E. Roda, B. Marchetti, F. Vincenzi, G. Zauli, P. A. Borea, C. A. Locatelli, K. Varani, M. Marti, 2023, "MAM-2201 acute administration impairs motor, sensorimotor, prepulse inhibition, and memory functions in mice: a comparison with its analogue AM-2201", Psychopharmacology

David L. Braff, Mark A. Geyer, Neal R. Swerdlow, 2001, "Human studies of prepulse inhibition of startle: normal subjects, patient groups, and pharmacological studies", Psychopharmacology

M. Marti, G. Talani, C. Miliano, S. Bilel, F. Biggio, J. Bratzu, M. Diana, M. A. De Luca, L. Fattore, 2021 "New insights into methoxetamine mechanisms of action: Focus on serotonergic 5-HT2 receptors in pharmacological and behavioral effects in the rat", Experimental Neurology

S. Bilel, M. Tirri, R. Arfè, S. Stopponi, L. Soverchia, R. Ciccocioppo, P. Frisoni, S. Strano-Rossi, C. Miliano, F. De-Giorgio2, G. Serpelloni, 2019, "Pharmacological and Behavioral Effects of the Synthetic Cannabinoid AKB48 in Rats", Frontiers in Neuroscience

M. Tirri, S. Bilel, R. Arfè, G. Corli, B. Marchetti, T. Bernardi, F. Boccuto, G. Serpelloni, F. Botrè, F. De-Giorgio, K. Golembiowska, M. Marti, 2022, "Effect of -NBOMe Compounds on Sensorimotor, Motor, and Prepulse Inhibition Responses in Mice in Comparison With the 2C Analogs and Lysergic Acid Diethylamide: From Preclinical Evidence to Forensic Implication in Driving Under the Influence of Drugs", Frontiers in Psychiatry

M. Tirri, L. Ponzoni, S. Bilel, Raffaella Arfè, D. Braida, M. Sala, M. Marti, 2020, "Acute DOB and PMA Administration Impairs Motor and Sensorimotor Responses in Mice and Causes Hallucinogenic E ects in Adult Zebrafish", Brain Sciences

STARTLE/PPI SYSTEM



The Ugo Basile Startle/PPI system is composed of:

- Isolation Cubicle (inside dimensions 39(d)x47(w)x49(h)cm) 46000-596, each with: dual visible/I.R. light 46000-325, speaker 46000-165, controller on-board 46000-105 and Stimulating/Recording Platform for mouse or rat (48006 or 48002). A 4-channel startle Box 48163-150 also included.
- The System needs to be completed with the included mouse or rat holder (select one when ordering).

Specifications - Ope	ration					
Sound Range	From 65 dB to 120 dE	3				
Pulse/pre pulse sound rar	nge From 65 dB to 120 dE	From 65 dB to 120 dB				
Sound Frequency Range	From 100 Hz to 18 KH	From 100 Hz to 18 KHz				
IR Light	0-100% (Environment	0-100% (Environment)				
Flash Light	0-100% (20K Lux) S/R	0-100% (20K Lux) S/R Platform)				
Data Export	.csv, .xlsx format	.csv, .xlsx format				
Specifications - Gen	eral					
Power Requirements	Universal input 100-2	Universal input 100-240 VAC, 50-60Hz, 60W				
Background Noise (Ventila	ation) from 0.1 to 3mA, in 0.	from 0.1 to 3mA, in 0.1mA steps				
Operating environment	10°C to 40°C; 5% to 9	10°C to 40°C; 5% to 95% RH (non-condensing)				
Physical	O uter size	I nternal size				
Dimensions Cubicle	52(d)x62(w)x58(h)cm	39(d)x47(w)x49(h)cm				
Ordering informatio	ons					
48163	Startle/PPI System for Mice, 1 anim	System for Mice, 1 animal system				
48263	Startle/PPI System for Mice, 2 animal system					
48463	Startle/PPI System for Mice, 4 anim	nal system				
48162	Startle/PPI System for Rats, 1 anim	ystem for Rats, 1 animal system				
48262	Startle/PPI System for Rats, 2 animal system					
48462	Startle/PPI System for Rats, 4 animal system					
48003-350	Startle PPI Small cylindrical mouse holder (up to 25g)					
48003-340	Startle PPI Medium cylindrical mouse holder (up to 50g)					
48003-330		arge cylindrical mouse holder (up to 75g)				
48003-360	Startle PPI XXL cylindrical mouse h					
48003-370	Startle PPI Standard cylindrical rat holder (up to 250g)					
48003-380	Startle PPI XL cylindrical rat holder (up to 800g)					
Extra warranty (standard 12 months + 12 months with product registration)						
48163-UBC12	UB Care 12 Additional hardware warranty extension 12 months for Startle/PPI system 1 animal system					
48163-UBC24	UB Care 24 Additional hardware warranty extension 24 months for Startle/PPI system 1 animal system					
48263-UBC12	UB Care 12 Additional hardware warranty extension 12 months for Startle/PPI system 2 animal system					
48263-UBC24	UB Care 24 Additional hardware warranty extension 24 months for Startle/PPI system 2 animal system					
48463-UBC12	UB Care 12 Additional hardware warranty extension 12 months for Startle/PPI system 4 animal system					
48463-UBC24	UB Care 24 Additional hardware warra	ditional hardware warranty extension 24 months for Startle/PPI system 4 animal system				

ugobasile.com

more than 40,000 citations in the main bibliographic search engines.



Ugo Basile SRL Via Giuseppe Di Vittorio, 2 21036 Gemonio (VA) ITALY Tel. +39 0332 744574 Get a quote: sales@ugobasile.com



Partner area